

DIAL Measurement of NO₂ Concentration in a Stack Plume

Nobuo Sugimoto, Nobuo Takeuchi and Michio Okuda

National Institute for Environmental Studies
Yatabe, Tsukuba, Ibaraki 305 Japan

Differential absorption lidar (DIAL) has been recognized as a promising technique for range-resolved remote measurement of pollutant gas concentration.

A DIAL system for NO₂ measurements using a flashlamp-pumped dye laser has been developed at the National Institute for Environmental Studies (NIES). In this paper, we describe the outline of the system and measurement results of NO₂ concentration in a stack plume.

A schematic of the system is shown in Fig.1, and the system specifications are shown in Table 1. The system is computerized so that wavelength switching, laser firing and data acquisition are performed automatically. Acquired data are stored on magnetic tape, and processed to obtain NO₂ concentration profiles after measurement.

A stack plume from the Energy Center of NIES was measured from the roof of the main building about 120 m away. The NO₂ concentration profile obtained from 60 pulse-pair averages is shown in Fig.2. The range resolution is 120 m. The profile shows a peak in the vicinity of the stack, and the concentration at this peak is ~ 0.17 ppm with a probable error of ~ 0.05 ppm. This concentration is an average over 120 m including the stack plume. If NO₂ is assumed to exist only in the plume (~ 2 m in diameter), the NO₂ concentration in the plume is then estimated to be ~ 10 ppm, which is a reasonable value.

We are now investigating the possibility of constructing a two-wavelength lasing system.

Table 1. System Specifications

Laser type	Flashlamp-pumped dye laser
Wavelength tuning	Etalon
Wavelength	λ_{on} 4631 Å λ_{off} 4658 Å
Dye	C1H
Pulse energy	~ 100 mJ
Pulse duration	600 ns
Repetition	0.5 pps
Telescope	0.15 m ϕ Newtonian
Transient recorder	DM 901 (Iwatsu) 10 ns minimum, 8 bits, 1024 segments
Minicomputer system	NOVA 02




Fig. 1 Schematic diagram of the System




Fig. 2 NO₂ concentration profile obtained from 60 pulse-pair averages